Änderungsnummer / Change Master: C2365_210

Objektverknüpfungen / Object Links:

JBX 40037492 000 01 MDR Co-Plus

Dokumentenstückliste / Document Structure:

Status		Responsible	Date	
IE	in Erstellung	DEHNINA	15.02.2021	
	in Erstellung			
AF	Freigabeanford.	DEHOLHS	15.02.2021	
FR	freigegeben	DEHNIPL	15.02.2021	

Product Data Sheet

1.0 N	ame of the product	Co-Plus [®]		
2.0 Pi	0 Product description			
2.1 D	escription	Co-Plus [®] is an elastic cohesive bandage made of polyamide, polypropylene and polyurethane coated with natural rubber latex.		
2.2 In	ntended purpose	Co-Plus [®] is intended for fixation of wound dressings of all types and sizes as well as devices such as splints, tubes and paddings, for short term application of compression and to provide general support. Co-Plus [®] is non-sterile, for single use on intact skin only. Intended users are healthcare professionals and patients. Use of the device is not restricted to specific populations.		
2.3 In	nstructions for use	Instructions for use not necessary for class I/Is/IIa products. Product is well-known and self-explanatory. Pictograms and symbols on folding box are guiding to a safe use.		
	/arnings and recautions for use	Caution: This product contains natural rubber latex which may cause allergic reactions. Do not re-use. Upper limit of temperature 30°C.		
2.5 C	ontraindications	Do not use Co-Plus [®] on patients with known intolerance to latex.		
2.6 Tı	ransport Precautions	Not applicable		
2.7 D	uration of application	The bandages can be used from one day to several days depending on the treatment needs but cumulative use can be longer than 30 days. For single use on intact skin only.		
re	nformation on eusable/ washable roducts	Do not re-use		
2.9 C	omposition	62.5% Natural rubber compound 35.3% Polyamide/Polypropylene 2.2% Polyurethane yarn		

essity

BSN documentation system JBM 40047786 000 02 PDS MDR Co-Plus

Released:

15.02.21-...

Seite 3 von 5

	Size	Color	Per flo w	Per folding box	Per shipper	Product code
	10cm x 3m (6.3m)	Flesh	1		18	72100-04
	15cm x 3m (6.3m)	Flesh	1		12	72100-05
	2.5cm x 3m (6.3m)	Flesh	5		60	72100-06
	5cm x 3m (6.3m)	Flesh	1		36	72100-08
	7.5cm x 3m (6.3m)	Flesh	1		24	72100-09
	10cm x 3m (6.3m)	White	1		18	72101-05
	5cm x 3m (6.3m)	White	1		36	72101-09
	7.5cm x 3m (6.3m)	White	1		24	72101-10
	10cm x 4.5m	White	1		16	72101-11 REPACK: 72101-14
	7.5cm x 4.5m	White	1		24	72101-12 REPACK: 72101-15
	5cm x 4.5m	White	1		36	72101-13 REPACK: 72101-16
	7cm x 3m	White	-	1	24	72101-27
2.10 Product	7cm x 3m	Flesh	-	1	24	72101-28
range	10cm x 3.5m	White	-	1	18	72101-29
	10cm x 3.5m	Flesh	-	1	18	72101-30
	10cm x 3m (6.3m)	3x Purple 3x Blue 3x Green 3x Pink 3x Red 3x Light Blue	1		18	72102-02
	7.5cm x 3m (6.3m)	4x Purple 4x Blue 4x Green 4x Pink 4x Red 4x Light Blue	1		24	72102-03
	5cm x 3m (6.3m)	6x Purple 6x Blue 6x Green 6x Pink 6x Red 6x Light Blue	1		36	72102-04
	10cm x 3m (6.3m)	Green	1		18	72102-05
	10cm x 3m (6.3m)	Blue	1		18	72102-23

BSN documentation system JBM 40047786 000 02 PDS MDR Co-Plus

2.11	Sizing Chart	Please see 2.10			
2.12	Storage conditions	Dry, clean and without exposure to direct sunlight. Max. storage condition: 30°C			
2.13	Shelf life/ Storage time	The product has a shelf life of 3 years, if stored appropriately			
3.0 C	hemical substances of s	of special concern to Essity			
A refer		stances that are of special concern and ices can be found in Annex A2 via the fo omponents.			
		Raw materials used in product formula	tion		
		Substance	Included in formulation		
		natural rubber Latex	Yes		
		Lanolin and its derivatives	no		
		Colophony	not intentionally added only rosin esters used		
		Colophony derivatives	not intentionally added only rosin esters used		
		Bisphenol A (BPA)	no		
		Polyvinylchloride (PVC)	no		
	Components	Methylmethacrylate	no		
		Butylacrylate	no		
		Microplastics (<5 mm)	no		
		Nanoparticles	no		
		Fragrance/parfum	no		
3.1 C		Antibiotics	no		
		Triclosan	no		
		Chlorhexidine	no		
		Polyhexanide	no		
		Raw materials used in the packaging			
		Substance	Included in formulation		
		Polyvinylchloride (PVC)	no		
		Natural rubber latex	no		
		Recycled material	no		
		Raw materials used in the product formulation and in the packaging may contain small amounts of the substances listed above and not marked as included in formulation as well as amounts of other substances not listed above. As trace level analysis is not performed for each batch BSN medical GmbH cannot provide proof of absence.			

BSN documentation system JBM 40047786 000 02 PDS MDR Co-Plus

History

Version/Date	Page /Item	Description of Change
01 / 24.06.2011	all	Set up new document 40014134 000 01
01 / 02.2021	all	Set up new document according to new template 40014134 000 01 will be replaced by 40047786
02 / 02.2021	Background data	Linked JBX